

# EVE

by Cryoskin

Inspired by K-Beauty Science  
Unlock the Power of Multi-Tech



CE FC FDA  
Listed

CRYOSKIN®

# Hydra-dermabrasion

## Beginning of skin care and build skin immunity

### Hydra-dermabrasion is the first step in skin care.

No matter how great the skin care products that you've applied, acne and inflammations will appear on the skin. The reason for this is the **Demodex folliculorum** and **Body waste**.

Demodex folliculorum is an enemy of the skin that destroys the skin balance by destroying collagen as well as producing acne and causing inflammation.

Leaving Demodex folliculorum and body waste in your skin is a shortcut to destroying the underlying skin balance.

Hydra-dermabrasion removes skin wastes and at the same time penetrates the skin with a solution suitable for all skin types.

Smoothens the skin texture inside and out, makes you experience a greater effect when undergoing other skin care procedures.



### How EVE works

As the exclusive solution "X-Stream" comes out from the tip, it rotates with a string radial principle, and by sucking waste materials, it is possible to remove wastes from the pores without scratching the skin.



#### Step1 Cleansing

Removing demodex folliculorum and skin wastes.



#### Step2 Extractions

Painless, automated extractions use vacuum suction to draw out impurities.



#### Step3 Peeling

AHA & PHA & LHA & Enzyme help the skin loosen dirt and debris from pores without irritation.



#### Step4 Hydration

Solution come out from the tips to nourish and protect the skin.

## Feature

#### Solution efficiency

Increased solution efficiency with the specially developed tips that minimizes the disadvantages of the existing vortex tips. It doesn't break down easily.

#### Powerful suction & adjustable levels

Minimize skin irritation and make the inside of the skin clean with a strong vacuum and penetrating power. Adjust the level of the vacuum power, suitable to users.

#### Ergonomic hand piece

Developed with an ergonomic design, it enables a more comfortable treatment.

#### Ease of use

- LCD touch screen
- LED GUIDE on tips
- Easy to plug and unplug solutions and sediment bottles

#### Available for all Skin types


Exclusive solution for Eve that contains LHA (4th generation), a newly developed ingredient from AHA, BHA, and PHA, and can be used for all skin types.

## What you can expect


- Cleansing
- Peeling
- Hydration
- Brightening
- Exfoliation (blackhead, whitehead, demodex folliculorum)
- Toning
- Acne
- Sebum control

## Various Tips


Large  
Small  
Cleaning  
Diamond




Small



Large



Scalp



Cleaning



Diamond




# Microcurrent / EMS

## The beginning and end of Anti-aging

Micro-current therapy is where Dr. Erwin Neher and Dr. Bert Sakmann won the 1991 Nobel Prize in Physiology and proved that there is a possibility of recovery when microcurrent is applied to damaged tissue. EVE has developed a microcurrent technique based on this theory. We pursue fundamental skin beauty.

### How EVE works

Microcurrent electrotherapy technology increases the production of adenosine triphosphate (ATP) and protein, which plays an important role in cellular energy metabolism, using a low current of  $\mu A$ , which is hardly felt by the human body, thereby promoting skin cell metabolism and lymphatic massage. Pain relief and skin elasticity induction technology.



### Feature

- **Non-invasive, painless** non-surgical way to reduce wrinkles and tighten skin in a safe, trustworthy and effective way.
- **Quick & Effortless** easy to learn, simple to operate, fast in delivering results.
- **NO down time**
- **Immediate & Long lasting effects** Eve provides fast and continuous anti-aging results.
- **Variations of Pre-set protocol programs** Easy to use even for beginners as the program is set for the purpose of using the device for skin problems such as anti-aging (by age), acne, etc.

### What you can expect

- Lifting
- Enhance definition
- Tighten and firm skin
- Stimulate collagen growth
- Drained out lymphatics
- Reduces puffiness
- Improves muscle tone
- Increases cellular activity
- Improved oxygenation of tissue and skin
- Pain relief
- Diminished pigmentation

### Various Hand piece

#### Ball type

The micro size of the ball intensifies the current and increases effectiveness.

#### T type

To lift broader areas such as the cheekbones, periorbital wrinkle and jaw line.

#### Comb Type


To energize hair follicles, promote blood circulation, diminish hair loss, reduce dryness and dandruff

#### Gloves

Ideal way to provide a more detail oriented microcurrent treatment

#### Electrode Patch (Option)

Easy, hands-free treatment targeting the fine lines and wrinkles



### • Storage Space



### • Solution

- New ingredient LHA & Enzyme included

### • Sensor

- Preventing overflow of sediment bottle
- Preventing malfunction

### • Easy to Plug and unplug

- No more wasting time to plug solution and sediment bottle.



### • Detachable hose

- Easy maintenance
- Preventing leakage and break down



## • Various Applicators



## • LED GUIDE

- BLUE LED LIGHTS (UV) : Targets and disrupts propionibacterium acnes, improves oily and congested skin
- WHITE LED LIGHTS : Guidelines for users



## • GUI (optimizing your treatment to the fullest)

- The advanced user interface is not only user friendly, but easy to adopt



## • EASY Pre-set protocol programs

- Total of 10 different pre-set protocols that allow you to easily provide customized treatment based on specific skin conditions or concerns in a controlled time setting.



# X-Stream

## AHA (Alpha Hydroxy Acid) Lactic Acid

AHA is a famous cosmetic ingredient derived from fruits, nuts, milk or sugars. AHA is water soluble and acts on the skin's surface. It acts as an exfoliator and retains water in the skin. It is useful for people with dry, rough skin such as in keratosis pilaris.



### AHA Feature

- Peeling thickest layer of dead skin cells
- Water soluble
- Inhibit pigment
- Hydrating
- Boost Collagen synthesis
- Antimicrobial
- Adjust PH

## PHA (Poly Hydroxy Acid) Lactobionic Acid

PHA is useful for people with sensitive skin such as rosacea or eczema because they have larger molecular structures. PHA can provide visible improvements in sun-damaged skin.

### PHA Feature

- Good for sensitive skin
- Soft exfoliation
- Great hydrating
- Anti-inflammatory
- Antioxidant
- Boost collagen synthesis
- Skin tone improvement
- Wrinkle improvement
- UV protection

## Enzyme Bromelain, Hyaluronic Acid, Niacinamide

Together, these enzymes deliver 8 major skin benefits at once: strengthens skin immune system, antioxidant, improves skin texture, promotes skin clarity, revitalizes dull complexion, hydrates, soothes and restores the natural skin's water and oil balance.



## LHA (Lipo Hydroxy Acid) Capryloyl Salicylic Acid

LHA is advanced version of BHA, one of the skincare industry's latest attempts at cracking the acne cure.

The major benefit of LHA, lies in its gentle, non-irritating mode of exfoliation.

LHA may also augment the anti-aging effects of tretinoin, an active form of vitamin A, by stimulating skin cell production.

### LHA Feature

- Best acne care treatment
- Gentle exfoliation
- Non-irritating
- Hydration
- Anti-aging
- Reduce redness
- Tight pores

# Pre-set protocol programs

EVE provides a skin enhancement program just for you.

Enjoy an explicit anti-aging effect with the EVE exclusive programs

## Multi-Step Programs

- Anti-Aging group by age
- Facial Acne/ seborrhea
- Rosacea
- Restoration post-destruction

## EMS Programs

- Complete comprehensive anti-aging care through EMS Special session.
- Stimulate and tone muscle to tighten the skin.
- Focus on deep skin areas where skin elasticity needs to be improved.

## Single-step programs

- **Balancing Skin Metabolism**

Balancing the mitochondrial membrane, thus improving skin cell metabolism

- **Alleviating Pain**


Alleviating pain from facial myalgias by electrical nerve stimulation

- **Lymphatic Drainage**

Improving blood and lymph flow  
Removing waste from the skin cells

- **Myo-stimulation**

Stimulating facial muscles (lifting of sagging skin)/ relaxing muscles in hyper tonus



# EVE

by Cryoskin

## Specification

<b>Size</b> (Net)	20 x 40 x 30 cm (Device) 36 x 40 x 105 cm (Cart)
<b>Size</b> (Packaging)	35 x 55 x 50 cm (Device) 48 x 18 x 76 cm (Trolley)
<b>Weight</b> (Packaging)	13 kg (Device) 8.5 kg (Cart)
<b>Suction power</b> (Hydra)	650mmHg
<b>Current</b> (Microcurrent)	2~500 UA
<b>Power</b> <b>Consumption</b>	100~240V, 50/60Hz, Max 100VA (24V/ 5A)

Cryoskin International Paris Headquarter

[contact@cryoskinparis.com](mailto:contact@cryoskinparis.com) | [www.cryoskin-international.com](http://www.cryoskin-international.com)